


Manufacturer
Steril

Model
Gemini



Vertical Laminar Flow Cabinet



Technical specifications:

- Laminar Airflow Cabinets are used for work with low-risk substances and materials and for all the situations where is requested to protect the product from dangerous effects due to uncontrolled diffusion of air-transported contaminants during its handling.
- Are designed in order to obtain a separated working area with a superior level of sterility .

Applications and uses:

- Normally used in laboratory of microbiology, molecular biology and, furthermore, used in innovative fields of scientific research.
- Since Vertical Laminar Flow Cabinets are not biological safety cabinets, they **CANNOT be used for pathogens handling.**